DENTSPLY
DETRAY

X e n o[®] III

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DENTSPLY
DETRAY

X e n o[®] III

Single Step Self-Etching Dental Adhesive

Xeno III¹ is a single step self-etching dental adhesive designed to bond light-cured restoratives to enamel and dentine.

Xeno III has the following features:

- The single step bonding system treats both enamel and dentin simultaneously with acid-conditioning, priming and bonding functions.
- No water rinsing after application.

COMPOSITION

Liquid A

2-hydroxyethyl methacrylate (HEMA)

Purified water

Ethanol

Butylated hydroxy toluene (BHT)

Highly dispersed silicon dioxide

Liquid B

Phosphoric acid modified methacrylate (Pyro-EMA)

Mono fluoro phosphazene modified methacrylate (PEM-F)

Urethane dimethacrylate

Butylated hydroxy toluene (BHT)

Camphorquinone

Ethyl-4-dimethylaminobenzoate

¹ Xeno III is the third generation of a self-etching dental adhesive developed by DENTSPLY Sankin and marketed in Japan under the brand Xeno CF II Bond.

INDICATIONS

Single step self-etching adhesive for **direct, light-cured** restorative materials.

CONTRAINDICATIONS

- Use with patients who have a history of severe allergic reactions to methacrylate resins or any other of the components.
- Direct or indirect pulpcapping.
- Use with chemical or dual-cure materials.

WARNINGS

1. Eye contact: Xeno III contains methacrylates which may be irritating to eyes. Before using this product wear protective glasses as well as covering the patient's eyes to protect from splash-

ing material. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

- 2. Skin contact:** Xeno III contains polymerisable monomers which can cause skin sensitisation (allergic contact dermatitis) in susceptible individuals. If contact with skin, immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water after contact. If skin rash and sensitisation or other allergic reaction occurs, discontinue use and consult a physician.
- 3. Oral mucosa contact:** Avoid contact with oral soft tissues. If accidental contact occurs, remove adhesive with cotton pellet and rinse as soon as the clinical situation allows this. After prolonged contact with the adhesive the alveolar mucosa may turn white and/or numb due to the contact of the adhesive. This is a temporary phenomenon and should disappear within a day. If it lasts longer, the patient should contact his dental professional.

4. **Ingestion:** Xeno III is acidic. Use caution from the patient accidentally swallowing the product. If accidental swallowing occurs, drink lots of water.

PRECAUTIONS

1. Xeno III is a light-cured material. Proceed immediately once the liquids have been placed into the mixing well or protect from ambient light.
2. Xeno III is less effective on unground, unprepared, or uncut enamel. Overfilling of restorative onto unetched, unground, or uncut enamel could cause marginal discolouration.

3. Avoid excessive contact of the adhesive with gingival tissues, as this may cause a transient whitening and/or numbness of the alveolar mucosa. (See Warnings section).

4. Interactions with dental materials

- 4.1 Insufficient data exist to support the use of Xeno III with indirect restorations.
- 4.2 Eugenol containing dental materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
- 4.3 If H₂O₂ has been used to clean the cavity, proper rinsing is essential. Higher concentration H₂O₂ may interfere with the setting of polymerisable material and should not be used prior to the application of Xeno III.

5. Xeno III bottles should be tightly closed immediately after use. Do not interchange caps.
6. This product is intended to be used only as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.

ADVERSE REACTIONS

1. Corneal damage may result with prolonged eye contact and/or exposure to Xeno III.
2. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals to Xeno III. (See Contra-indications, Warnings, and Precautions sections).

3. Transient whitening and/or numbness of the alveolar mucosa may occur after contact with Xeno III.

STEP-BY-STEP INSTRUCTIONS

1. Cleaning and Preparation

Cavity cleanliness is paramount for the development of adhesion.

In cases where no cavity preparation has been made, clean the tooth surface with a rubber cup and pumice or a prophylaxis paste (e.g. Nupro®). Preparing a fresh surface with a finishing bur will significantly increase bond strength to enamel. In all situations where maximum retention is required instrument enamel (e.g. with a finishing diamond bur) to achieve optimal bonding results.

Wash surface thoroughly with air/water spray.

Remove rinsing water by blowing gently with an air syringe or blot-dry with a cotton pellet.

Do not desiccate the dentine structure.

2. Pulp Protection

For direct or indirect pulp capping cover the dentine close to the pulp (< 1 mm) with a hard-setting calcium hydroxide liner (e.g. Dycal[®]), leaving the rest of the cavity surface free for bonding with Xeno III.

3. Application

3.1 Dispense an equal amount of Liquid A and Liquid B into a clean mixing well.

For optimum dosage, proceed as follows:

Liquid A

Shake bottle two to three times prior to use.

Incline bottle. If liquid does not drop, squeeze.

After use, close bottle with green cap.

Liquid B

Incline bottle; wait for approx. two seconds. Squeeze until liquid drops.

After use, close bottle with black cap.

3.2 Mix liquid thoroughly for approx. five seconds with the applicator tip supplied.

If protected from ambient light, the mixed material stays functional for up to 60 minutes!

3.3 With the applicator tip provided, apply generous amounts of Xeno III to thoroughly wet all cavity surfaces.

- 3.4 Leave undisturbed for at least 20 seconds.
- 3.5 Uniformly spread the adhesive using a gentle stream of oil free air for at least 2 seconds until there is no more flow of the adhesive to ensure proper removal of solvent. Avoid thinning the adhesive layer by excessive air-drying.
Do NOT water rinse.
- 3.6 Cure² the adhesive with a light-curing unit for at least 10 seconds.
- 3.7 Place restorative immediately following manufacturer's instructions.

² Check curing light for minimum curing output of at least 500 mW/cm².

STORAGE

Keep out of direct sunlight and protect from moisture.
To be stored at temperatures between 2 °C and 24 °C.
Inadequate storage conditions will shorten the shelf life and may lead to a malfunction of the product.
Close caps tightly after use.

BATCH CODE LOT AND EXPIRATION DATE

The batch code should be quoted in all correspondence, which requires identification of the product.
Do not use after expiration date.

If you have any questions, please contact:

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© DENTSPLY DeTrey 2003-05-08
DENTSPLY DeTrey GmbH