---

**DENTSPLY**  
DETRAY

*Dyract<sup>®</sup>flow*

**Flowable Compomer Restorative**

<b>Directions for Use</b> _____	<b>2</b>	<b>Οδηγίες χρήσης</b> _____	<b>72</b>
<b>Gebrauchsanweisung</b> _____	<b>12</b>	<b>Bruksanvisning</b> _____	<b>82</b>
<b>Mode d'emploi</b> _____	<b>22</b>	<b>Brugsanvisning</b> _____	<b>92</b>
<b>Gebbruiksaanwijzing</b> _____	<b>32</b>	<b>Инструкция по применению</b> _____	<b>102</b>
<b>Istruzioni per l'uso</b> _____	<b>42</b>	<b>Instrukcja użytkowania</b> _____	<b>112</b>
<b>Instrucciones de uso</b> _____	<b>52</b>	<b>Návod k použití</b> _____	<b>122</b>
<b>Instruções de utilização</b> _____	<b>62</b>	<b>Használati utasítás</b> _____	<b>132</b>

**DENTSPLY**  
DETRAY

*Dyract<sup>®</sup>flow*

**Flowable Compomer Restorative**

**Dyract® flow** is a universal compomer restorative material with *flow* characteristics which make it ideal for small cavities in anterior and posterior teeth.

**Dyract flow** adapts to the cavity walls without the use of hand instruments.

**Dyract flow** is delivered in small syringes for direct intra-oral application.

**Dyract flow** is available in 5 Vita<sup>®1</sup> shades in a regular opacity and additionally in a translucent and an opaque shade.

**Dyract flow** is used following application of Prime&Bond<sup>®</sup> NT, a universal self-priming dental adhesive designed to bond the restorative to enamel and dentine.

**Caution:** For dental use only.

## COMPOSITION

### Dyract flow

- Strontium-alumino-fluoro-silicate glass
- Highly dispersed silicon dioxide
- Ammonium salt of PENTA (dipentaerythritol penta acrylate monophosphate) and N,N-dimethyl aminoethyl methacrylate
- Carboxylic acid modified methacrylate macromonomers
- Diethylene glycol dimethacrylate (DGDMA)
- Camphorquinone
- Ethyl-4-dimethylaminobenzoate
- 2-Hydroxymethoxybenzophenone
- Butylated hydroxy toluene (BHT) and other stabilisers
- Iron pigments
- Titanium dioxide

---

<sup>1</sup> Vita is a registered trademark of Vita Zahnfabrik H. Rauter GmbH & Co.

### **Prime&Bond NT**

- Di- and trimethacrylate resins
- Urethandimethacrylate (UDMA)
- Highly dispersed silicon dioxide
- PENTA (dipentaerythritol penta acrylate monophosphate)
- Camphorquinone
- Ethyl-4-dimethylaminobenzoate
- Butylated hydroxy toluene (BHT)
- Cetylamine hydrofluoride
- Acetone

### **INDICATIONS**

- Minimal invasive restorations – this comprises:
  - Small, non occlusal stress-bearing restorations in anterior and posterior teeth
  - Extended pit & fissure sealings
  - Shallow class V

- Margin repair
- Cavity lining
- Adhesive cementation of indirect, light-transmitting porcelain or composite inlays and onlays.

### **CONTRAINDICATIONS**

#### **Prime&Bond NT**

- Direct or indirect pulp capping
- Use on patients with a known allergy to dimethacrylate resins or any other of the components.

#### **Dyract flow**

- Use on patients with a known allergy to dimethacrylate resins or any other of the components
- Permanent restorations of occlusal stress-bearing cavities in permanent posterior teeth

- When moisture control cannot be guaranteed during application (deep subgingival cavity margins)
- Core build-up for full ceramic crowns
- Luting of crowns or veneers
- Luting of restorations that are impenetrable by dental curing light or the walls of which exceed a thickness of 5 mm.

## WARNINGS

1. Prime&Bond NT and Dyract *flow* contain methacrylates that may be irritating to the eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
2. Prime&Bond NT and Dyract *flow* may cause sensitisation by contact with skin or mucous membranes in susceptible persons. After accidental contact, wash immediately with plenty of soap and water or rinse with plenty of water.  
If sensitisation has occurred, discontinue use.

3. Prime&Bond NT contains acetone. Acetone is highly flammable. Keep away from sources of ignition – no smoking. Do not breathe vapour. Take precautionary measures against static discharges.

## PRECAUTIONS

Avoid Prime&Bond NT saturating gingival retraction cord. If Prime&Bond NT soaks into the cord, it may set hard and bond the cord to the underlying tooth surface making removal difficult.

## INTERACTIONS WITH DENTAL MATERIALS

Eugenol containing dental materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.

If H<sub>2</sub>O<sub>2</sub> has been used to clean the cavity, proper rinsing is essential. Higher concentration H<sub>2</sub>O<sub>2</sub> may interfere with the setting of polymerisable material and should not be used prior to the application of Prime&Bond NT.

Prolonged and intensive contact with acetone-containing products may lead to minute dissolution of the outermost surface of calcium hydroxide materials. This has no detrimental effect on the adhesion to the cavity walls.

## **ADVERSE REACTIONS**

The following adverse reaction has been associated with the use of acetone solutions and acrylate monomers:

- Reversible inflammatory changes of the oral mucosa after accidental contact.

## **STEP-BY-STEP INSTRUCTIONS**

### **I. DIRECT RESTORATIONS**

#### **1. Shade Selection and Cleaning**

Shade selection should be made prior to the restorative procedure whilst the teeth are hydrated. Remove any extraneous plaque or surface stain

with a rubber cup and pumice or a prophylaxis paste like Nupro®. Use the Dyract *flow* shade guide which contains samples of original Dyract *flow* restorative. The colour coding dot on the shade guide matches the coloured ring on the label of the syringe.

Alternatively, a Vita Lumin® Vacuum shade guide may be used. The Dyract *flow* shade corresponds to the central part of the respective Vita tooth.

#### **2. Cavity Preparation**

In all classes of cavity this may be kept to the minimum required for caries removal.

In cases where no cavity preparation has been made refresh the surface with a finishing bur.

Wash surface thoroughly with air/water spray.

Remove rinsing water by blowing gently with an air syringe or blot-dry with a cotton pellet.

Do not desiccate the dentine structure.

### 3. Pulp Protection

For direct or indirect pulp-capping cover the dentine close to the pulp with a hard-setting calcium hydroxide liner (e.g. Dycal®) leaving the rest of the cavity surface free for bonding with Prime&Bond NT.

### 4. Conditioning of dentine and enamel

For most restorative procedures with Dyract *flow* it is not necessary to condition the prepared tooth. In this case, **proceed to step 5.**

If Dyract *flow* is used as a cavity liner under a restorative material, follow the recommended conditioning procedure for the respective restorative, e.g. DeTrey® Conditioner 36 for Esthet•X® and other composites.

If surfaces are to be treated which are placed exclusively in enamel (e.g. use as extended pit and fissure sealant, preventive resin restoration, or for approximal slot preparations) and for cavities with bevelled enamel

margins, an acid conditioning procedure should be carried out. Follow the directions for use for the respective conditioner.

Once the surfaces have been properly conditioned, they must be kept uncontaminated. If salivary contamination occurs, thoroughly clean with forceful water spray and repeat the acid conditioning procedure.

### 5. Application of Prime&Bond NT

5.1 Dispense Prime&Bond NT directly onto a fresh Applicator Tip<sup>2</sup> or onto a disposable brush. Alternatively, dispense into a fresh DENTSPLY Applicator Dish<sup>2</sup> or standard dappen dish.

5.2 Immediately apply ample amounts of Prime&Bond NT to thoroughly wet all cavity surfaces.

These surfaces should remain fully wet for 20 seconds and may necessitate additional application of adhesive.

---

<sup>2</sup> DENTSPLY Applicator Dish and Applicator Tips are available from your dental dealer.

Remove solvent by blowing gently with air from a dental syringe for at least 5 seconds. Surface should have a uniform, glossy appearance. If not, repeat application and air-dry.

5.3 Light-cure for a minimum of 10 seconds<sup>3</sup>. Ensure uniform exposure of all cavity surfaces.

5.4 Immediately place Dyract *flow* over the cured Prime&Bond NT.

## 6. Placement of Dyract *flow*

6.1 Remove cap from the end of the syringe. To ensure free flow of material from syringe, express a small amount onto a pad.

6.2 Attach disposable needle tip to end of the syringe. Turn tip clockwise 1/4 to 1/2 turn to assure that it is fully seated. Tug on tip to be sure that it is locked into the collar of the syringe.

6.3 Dyract *flow* should flow freely with gentle pressure. DO NOT USE EXCESSIVE FORCE. If more than gentle pressure is required, remove from patient field and check for obstruction.

6.4 Dispense Dyract *flow* directly into the cavity preparation. In deep cavities, incremental placement and curing (in 2 mm layers or less) is recommended to minimise polymerisation shrinkage.

6.5 DISCARD NEEDLE IMMEDIATELY AFTER USE. CLOSE SYRINGE WITH ORIGINAL CAP<sup>4</sup>. Do not store syringe with dispensing tip in place. STORE ONLY WITH ORIGINAL CAP.

6.6 It is recommended to pull back slightly on the syringe plunger after use to prevent excessive flow of material.

## 7. Curing

Cure each increment separately with a VLC dental polymerisation unit for at least 20 seconds and 40 seconds for O-A3 and A4, respectively. The tip of the light guide should be held as close as possible to the restoration during curing.

---

<sup>3</sup> Check curing light for minimum curing output of at least 500 mW/cm<sup>2</sup>.

<sup>4</sup> Cap can be stored in the respective compartment of the starter pack.

**Important:** Be sure to expose each area of the entire restoration to the curing light. Additionally, the restoration should be cured through lingual or buccal enamel walls.

## **8. Finishing**

Finish immediately after curing. Gross excess material may be removed with fluted finishing burs or diamonds. Finishing is best achieved by using Enhance™ Finishing and Polishing Discs and interproximal finishing and polishing strips. A high final lustre can be obtained by applying Prisma® Gloss™ and Prisma Gloss Extrafine Polishing Pastes.

## **II. INDIRECT RESTORATIONS**

### **1. Pre-Treatment of Restoration**

The internal surfaces should be treated according to the manufacturers' recommendations.

### **2. Moisture Control**

Surface cleanliness is paramount for the development of adhesion.

### **3. Enamel and Dentine Pre-Treatment**

Clean enamel and dentine prepared during a previous appointment with a rubber cup and pumice or a cleaning paste such as Nupro Prophylaxis Pastes. Wash thoroughly with water spray and then air-dry. Clean freshly cut enamel and dentine with water spray and then air-dry.

### **4. Conditioning and Application of Prime&Bond NT**

For conditioning of enamel and dentine and application of Prime&Bond NT, follow the instructions for use described for direct restorations.

### **5. Dispensing of Dyract *flow***

Dispense Dyract *flow* onto the restoration and/or onto the tooth preparation.

## 6. Placement of Restoration

Place restoration. Apply pressure and then release to allow the hydraulic pressure to dissipate. Repeat this procedure at least one more time and ensure that the restoration is completely seated.

Remove as much excess as possible before the compomer is light-cured. If any excess compomer remains after curing, remove with curette, scalers or finishing burs/discs or strips.

## 7. Light-Curing

Expose all surfaces of the restorations to the dental curing light for 40 seconds each. Depending on the number of surfaces of the restoration, the diameter of the curing light outlet and the size of the tooth, this results in up to 5 curing cycles per restoration.

## 8. Occlusal Adjustments and Finishing

Make occlusal adjustments and finish the margins with finishing burs. Remove excess from proximal surfaces with diamond strips or with the diamond tips for the EVA contra-angle handpiece.

Use polishing discs (Enhance) on inlay margins where accessible. Use polishing pastes (Prisma Gloss and Prisma Gloss Extrafine) additionally or alternatively. Work interdentally with aluminium oxide strips.

## STORAGE

The Prime&Bond NT bottle and Dyract *flow* syringes should be tightly closed immediately after use.

Keep out of sunlight.

Keep Prime&Bond NT in a well ventilated place.

To be stored at temperatures between 10 °C and 24 °C.

Humidity can adversely affect the properties of unsealed Dyract *flow* syringes. Therefore keep syringes sealed in their blister pack until use.

Under normal ambient conditions, unsealed syringes stay usable for about 5 months.

Under high humidity (90% r.h.) unsealed syringes stay usable for at least 3 weeks.

**BATCH NUMBER LOT AND EXPIRY DATE **

The batch number should be quoted in all correspondence, which requires identification of the product.

Do not use after expiry date.

If you have any questions, please contact:

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© DENTSPLY DeTrey 2002-03-21
**DENTSPLY DeTrey GmbH**