DENTSPLY
DETREY

NRC™

Non-Rinse Conditioner

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Conditionneur sans rinçage

Conditioner zonder spoelen

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Grabado ácido sin lavado

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NRC™

Non-Rinse Conditioner

NRC™ is a self-priming conditioning liquid used in combination with Prime&Bond®NT Nano-Technology Dental Adhesive for Dyract® compomer restorations where enhanced adhesion to enamel is required.

NRC offers a fast way to condition and prime enamel and dentine in an easy application procedure without the need to rinse with water.

Caution: For dental use only.

COMPOSITION

Organic acids/monomers in aqueous solution.

INDICATIONS

NRC is to be used prior to application of Prime&Bond NT and Dyract compomer materials in the following situations:

- Stress-bearing class I and II permanent Dyract®AP restorations
- Class IV Dyract AP restorations
- Cavities where bevelled enamel margins are preferred for aesthetic reasons.

NRC may be used for all other Dyract applications, if desired.

CONTRAINDICATIONS

Application to dentine close to the pulp (less than 1 mm).
Use in combination with adhesives other than Prime&Bond NT.

WARNINGS

NRC contains organic acids which may cause burns. Avoid contact with oral tissues, eyes and skin. If accidental contact occurs, flush affected area with copious amounts of water. In case of contact with the eyes, immediately rinse with plenty of water and seek medical attention.

INTERACTIONS WITH DENTAL MATERIALS

Some liners and bases may be etched superficially by NRC. In general, this does not impair their barrier function.

STEP-BY-STEP INSTRUCTIONS

1. Cleaning

Cavity cleanliness is paramount for the development of adhesion.

In cases where no cavity preparation has been made, clean the tooth surface with a rubber cup and pumice or a prophylaxis-paste like Nupro®. Preparing a fresh surface with a finishing bur will significantly increase bond strength to enamel.

Wash surface thoroughly with air/water spray.

Remove rinsing water by blowing gently with an air syringe or blot-dry with a cotton pellet.

Do not desiccate the dentine structure.

2. Pulp Protection

For direct or indirect pulp-capping protect the dentine close to the pulp (< 1 mm) with a hard-setting calcium hydroxide liner (e.g. Dycal®) leaving the remaining cavity surface free for bonding with Prime&Bond NT.

3. Application of NRC

- 3.1 Dispense NRC into a DENTSPLY Applicator Dish¹ or standard dappen dish.
- 3.2 Apply sufficient amounts of NRC with an Applicator Tip¹ or disposable brush to enamel and dentine. Leave undisturbed for 20 seconds. **Do not rinse.**
- 3.3 Remove excess NRC by blowing gently with an air syringe or blot dry with a cotton pellet. Do not desiccate the dentine structure.

Once the surfaces have been properly treated, they must be kept uncontaminated. If salivary contamination occurs, thoroughly clean with forceful water-spray and repeat the application of NRC.

4. Application of Prime&Bond NT

One Layer of Prime&Bond NT is applied:

- 4.1 Dispense Prime&Bond NT directly onto a fresh Applicator Tip or onto a disposable brush. Alternatively, dispense into a fresh DENTSPLY Applicator Dish or standard dappen dish.
- 4.2 Immediately apply ample amounts of Prime&Bond NT to thoroughly wet all tooth surfaces. This surface should be saturated which may necessitate additional application of Prime&Bond NT.
- 4.3 Leave the surface undisturbed for 20 seconds.
- 4.4 Remove solvent by blowing with air from a dental syringe for at least 5 seconds. Surface should have a uniform, glossy appearance. If not, repeat steps 4.2 to 4.4.
- 4.5 Light-cure for a minimum of 10 seconds². Ensure uniform exposure of all cavity surfaces.
- 4.6 Immediately place Dyract compomer over the cured Prime&Bond NT.

STORAGE

Not to be stored at temperatures exceeding 24 °C.

BATCH NUMBER AND EXPIRY DATE

The batch number should be quoted in all correspondence which requires identification of the product.

Do not use after expiry date.

If you have any questions, please contact:

Manufacturer:
DENTSPLY DeTrey GmbH
De-Trey-Str. 1
78467 Konstanz
GERMANY
Phone +49 (0) 75 31 5 83-0

Distributor:
DENTSPLY Limited
Hamm Moor Lane
Addlestone, Weybridge
Surrey KT15 2SE
Phone (0 19 32) 85 34 22

© DENTSPLY DeTrey 1998-04

¹ DENTSPLY Applicator Dish and Applicator Tips are available from your dental dealer.

² When using a high-performance unit such as the curing-lights manufactured by DENTSPLY, a curing time of 10 seconds is sufficient. For curing lights with an output lower than 300 mW/cm², a curing time of 20 seconds should be observed.

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