**POGO™**

**One Step Diamond Micro-Polisher**

Polishing System  
Sistema de Pulido  
Instruments à polir  
Polier-System  
Sistema per lucidatura

Sistema de Polimento  
Polishing System  
Polijststysteem  
Polishing System  
Σύστημα στίλβωσης

---

**DENTSPLY**  
CAULK

**CE**  
0120

***Directions  
For Use***

<b>Directions for Use</b>	<b>2</b>
<b>Instrucciones de uso</b>	<b>6</b>
<b>Mode d'emploi</b>	<b>10</b>
<b>Gebrauchsanweisung</b>	<b>14</b>
<b>Istruzioni per l'uso</b>	<b>19</b>
<b>Instruções de utilização</b>	<b>23</b>
<b>Brugsanvisning</b>	<b>27</b>
<b>Gebruiksaanwijzing</b>	<b>31</b>
<b>Bruksanvisning</b>	<b>35</b>
<b>Οδηγίες χρήσης</b>	<b>39</b>

## **Directions for Use — English**

The PoGo™ Polishers (Discs, Cups and Points) are pre-mounted, single-use diamond impregnated cured urethane dimethacrylate resin polishing devices designed for use in the final polishing of all accessible composite resin restorations. When used as directed, the PoGo™ Polishing System will yield a high surface luster on hybrid composites such as Esthet•X® Micro Matrix, TPH® Spectrum®, Prisma® TPH®, and Prisma® AP.H® Composites, SureFil™ High Density Restorative Material, Quixfil™ and microfilled composite restoratives. The PoGo™ Polishing System is also effective in final polishing of compomer and flowable composite and compomer restoratives such as Esthet•X® Flow, X-Flow™, Dyract® AP, Dyract® eXtra and Dyract®flow.

Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

### **INDICATIONS**

PoGo™ Polishing Discs, Cups and Points are indicated as the final step in polishing all properly contoured and finished composite resin and compomer restorations.

### **COMPOSITION**

Polymerized Urethane Dimethacrylate Resin  
Fine diamond powder  
Silicon Oxide  
Plastic Latch-type mandrel

## CONTRAINDICATIONS

1. Do not reuse. PoGo™ Polishers are indicated as a single use only, disposable item. Attempts to sterilize, disinfect and/or reuse these items can adversely affect their performance.
2. PoGo™ POLISHER SHOULD NOT BE USED IN HIGH-SPEED HANDPIECES.

## WARNINGS

1. PoGo™ Polishing Disc, Cups and Points are designed for use without water. As with any rotary instrument, heat will build up with prolonged contact. Use with light intermittent pressure. Application of an air cooling stream during polishing assists with residue removal.
2. Use of the PoGo™ Polishers with excessive lateral pressure may dislodge polisher from the mandrel or cause quick deterioration and decreased effectiveness. This type of excess pressure may cause gouging of the restoration surfaces.
3. Care should be taken to have good control on any rotary instrument to protect patients from injury. Contact with soft tissues (gingiva, mucosa) may result in injury to the tissue. If contact occurs, thoroughly wash the affected area with water and seek appropriate medical attention.

## PRECAUTIONS

1. This product is intended to be used as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
2. Use in a well ventilated area. Avoid inhaling airborne particulates.
3. Protect adjacent restorations from accidental contact.
4. **Storage & Shelf Life: "Dry Storage is Recommended". Store the PoGo™ Polishers Discs between 65°F-80°F (18°C-25°C) and in relative humidity of 50% ± 20%. Avoid direct sunlight. Do not use after expiry date.**

## **ADVERSE REACTIONS**

Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.

## **STEP BY STEP INSTRUCTION FOR USE**

**STEP 1:** Place restoration or composite veneer as per standard techniques. Complete gross reduction of excess and general outline form of the restoration using carbide finishing burs (i.e. Prisma® Finishing Burs), Enhance® Finishers or diamond finishing burs.

**Technique Tip:** Prior to proceeding to the PoGo™ Polishing System (disc, cup, point) restoration surface should be finished to a desired contour. Surface should be smooth and defect free.

**STEP 2:** For best results in larger restorations, the use of an intermediate finishing instrument such as an Enhance® Finisher (disc, cup, point) in a conventional contra-angled handpiece will aid in removing scratches and establishing the desired smooth, defect free, matte surface. **NOTE: ENHANCE® FINISHERS MUST NOT BE USED IN HIGH SPEED HANDPIECES.** In smaller restorations, acceptable surface contours and suitable surface finish may be achieved by the appropriate contouring/finishing instrument. If so, proceed to Step 3.

**STEP 3:** Insert a PoGo™ Polishing Disc, Cup or Point into a conventional speed contra-angled handpiece. **NOTE: PoGo™ POLISHERS MUST NOT BE USED IN HIGH SPEED HANDPIECES.** Begin polishing by applying light intermittent pressure at moderate speed to the restoration surface. When using the disc, the broad surface of the PoGo™ Polisher should be used in contact with the restoration to optimize its effectiveness. To increase surface luster, decrease the pressure applied to the restoration surface using a light buffing motion.

Technique Tip:

1. PoGo™ Polishing Discs are ideal for the facial surfaces and established line angles.
2. PoGo™ Polishing Points are ideal for posterior occlusal and concave lingual surfaces.
3. PoGo™ Polishing Cups are well suited for cervical and proximal line angle areas.

#### **STERILIZATION AND DISINFECTION**

Do not attempt to sterilize and/or reuse the PoGo™ Polishing System (disc, cup or point.) Properly dispose after use.

#### **LOT NUMBER AND EXPIRATION DATE**

1. Do not use after expiration date. ISO standard used: “YYYY/MM”
2. The following numbers should be quoted in all correspondences.
  - Reorder number
  - Lot number on package
  - Expiration date on package

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©2002 DENTSPLY International Inc. All Rights Reserved.
Manufactured in the U.S.A. by:
DENTSPLY DeTrey GmbH
PRINTED IN U.S.A.