---

**DENTSPLY**  
DETREY

**Posterior Restorative**

## Standard

<b>Directions for Use</b> _____	<b>2</b>
<b>Gebrauchsanweisung</b> _____	<b>8</b>
<b>Mode d'emploi</b> _____	<b>14</b>
<b>Istruzioni per l'uso</b> _____	<b>20</b>
<b>Instrucciones de uso</b> _____	<b>26</b>

## Special

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**DENTSPLY**  
DETREY

**Posterior Restorative**

**QuixFil** is a resin-based restorative material specially designed for posterior use.

**QuixFil** allows a simplified and fast application technique without compromising safety and efficacy.

**QuixFil** utilises a new filler technology resulting in a unique high filler load as well as fluoride release.

**QuixFil** exhibits advanced performance regarding abrasion resistance and physical properties which are of special importance for occlusal stress-bearing restorations.

**QuixFil** is pre-dosed in Compules® tips for direct intra-oral application.

**QuixFil** is available in one universal shade.

**QuixFil** is used following application of Xeno® III Single Step Self-Etching Dental Adhesive, Prime&Bond® NT, or other adhesives designed for use with QuixFil.

**Note:** To be applied by a dentist only.

## **COMPOSITION**

### **QuixFil**

Urethane dimethacrylate (UDMA)

Triethyleneglycol dimethacrylate (TEGDMA)

Di- and trimethacrylate resins

Carboxylic acid modified dimethacrylate resin

Butylated hydroxy toluene (BHT)

UV stabiliser

Camphorquinone

Ethyl-4-dimethylaminobenzoate

Silanated strontium aluminum sodium fluoride phosphate silicate glass

The filler : matrix ratio is 85.5 : 14.5 by weight and 66.4 : 33.6 by volume, respectively.

## **INDICATIONS**

QuixFil is indicated for Class I and II cavities when safety and performance are to be combined with a fast and efficient application.

## **CONTRAINDICATIONS**

- Use in patients with a known allergy to methacrylate resins or any other of the components.
- When contamination with saliva, blood, etc. cannot be avoided.

## **WARNINGS**

1. QuixFil contains methacrylates that may be irritating to the eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

2. QuixFil may cause sensitisation by contact with skin or mucous membranes in susceptible persons. After accidental contact, wash immediately with plenty of soap and water or rinse with plenty of water. If sensitisation has occurred, discontinue use.

## **STEP-BY-STEP INSTRUCTIONS**

### **1. Preparation**

Clean tooth to be prepared with a rubber cap and a prophylaxis paste like Nupro®.

Preparation may be kept to the minimum required for caries removal. Finishing of cavity margins with a fine finishing bur will contribute to marginal quality.

After preparation, wash cavity thoroughly with air/water spray.

Remove rinsing water by blowing gently with an air syringe or blot-dry with a cotton pellet.

Do not desiccate the dentine structure.

## 2. Pulp Protection

For direct or indirect pulp-capping protect the dentine close to the pulp (< 1 mm) with a hard-setting calcium hydroxide liner (e.g. Dycal®), leaving the remaining cavity surface free for bonding with the adhesive.

## 3. Conditioning and application of adhesive

See directions for use for Xeno III Single Step Self-Etching Dental Adhesive.

## 4. Placement of QuixFil

Insert Compules tip into the notched opening of the applicator gun barrel.

Dispense QuixFil directly into the cavity preparation. Due to the high depth of cure and low shrinkage for QuixFil, placement in layers up to 4 mm thick is possible.

## 5. Curing

Cure each increment separately with a light curing unit:

- 20 seconds for units with an output of 500 to 800 mW/cm<sup>2</sup>.
- 10 seconds for units with an output of 800 mW/cm<sup>2</sup> and more.

The tip of the light guide should be held as close as possible to the restoration during curing.

**Important:** Be sure to expose each area of the entire restoration to the curing light.

Additionally, the restoration should be cured through lingual or buccal enamel walls.

## 6. Finishing

Begin finishing immediately after curing. Gross excess material may be removed with fluted finishing burs or diamonds. Finishing is best achieved by using Enhance™ Finishing and Polishing Discs and interproximal finishing and polishing strips. In patients with an adequate oral hygiene, the final high lustre of the restoration develops with use.

## **MAINTENANCE OF APPLICATOR GUN**

The applicator gun is sterilisable by autoclave or cold sterilisation solution following the manufacturers' instructions.

It is recommended that the applicator gun be disassembled for assured sterilisation. Partially close the applicator gun and place thumb under the rear portion of the hinge. Push upward and lift hinge separating the applicator gun, exposing the plunger. Remove residual restorative with a soft paper tissue and a suitable solvent (70% alcohol).

To reassemble, insert plunger into applicator gun barrel, press components together and snap hinge mechanism in place.

## **STORAGE**

QuixFil Compules tips should be tightly closed immediately after use. Keep out of sunlight.

To be stored at temperatures between 10 °C and 24 °C.

Humidity can adversely affect the properties of unsealed Compules tips. Therefore keep Compules tips sealed in their blister pack until use. Un-

der normal ambient conditions, unsealed Compules tips stay usable for about 4 weeks.

## **BATCH NUMBER LOT AND EXPIRY DATE**

The batch number should be quoted in all correspondence, which requires identification of the product.

Do not use after expiry date.

If you have any questions, please contact:

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**© DENTSPLY DeTrey 2002-04-29**
**DENTSPLY DeTrey GmbH**