---

**DENTSPLY**  
DETREY

SureFil™

**High Density Posterior Restorative**

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# SureFil™

## High Density Posterior Restorative

**SureFil™** High Density Posterior Restorative is a visible light activated, radiopaque restorative material designed for stress-bearing posterior restorations of primary and permanent teeth. It is to be used with Prime&Bond® NT Nano-Technology Dental Adhesive. The handling characteristics allow rapid, bulk placement and assist with establishment of properly located and shaped proximal contact areas. This one-component, Visible Light Cured SureFil is packaged in syringes and in individual dose, light opaque cups.

**Caution:** For dental use only.

### COMPOSITION

#### DeTrey® Conditioner 36:

- Phosphoric acid
- Highly dispersed silicon dioxide
- Detergent
- Pigment
- Water

#### Prime&Bond NT:

- Di- and trimethacrylate resins
- Functionalised amorphous silica
- PENTA (dipentaerythritol penta acrylate monophosphate)
- Photoinitiators
- Stabilisers
- Cetylamine hydrofluoride
- Acetone

#### SureFil:

- Inorganic fillers
- Resin matrix

### INDICATIONS

1. SureFil is indicated for cavity classes I, II, in posterior teeth.
2. SureFil may be used as a direct restorative material and for the fabrication of inlays and onlays.

### CONTRAINDICATIONS

#### Prime&Bond NT:

Direct or indirect pulpcapping.

#### Prime&Bond NT / SureFil:

SureFil and Prime&Bond NT Adhesive are contraindicated for use with patients that have a history of severe allergic reaction to methacrylate resins.

## WARNINGS

1. Prime&Bond NT contains methacrylates which may be irritating to skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of soap and water. The product may cause sensitisation by skin contact in susceptible persons. If skin sensitisation occurs discontinue use.
2. Avoid contact of Prime&Bond NT with mucous membranes. After accidental contact, wash and rinse with plenty of water.
3. Prime&Bond NT contains acetone. Acetone is highly flammable. Keep away from sources of ignition – no smoking. Do not breathe vapour. Take precautionary measures against static discharges.
4. SureFil contains polymerizable monomers which can cause skin sensitization (allergic contact dermatitis) in susceptible individuals. Wash thoroughly with soap and water after contact. If skin sensitization or other allergic reaction occurs, discontinue use.
5. DeTrey Conditioner 36 contains 36% phosphoric acid. Causes burns. Avoid contact with oral tissues, eyes and skin. If accidental contact occurs, flush affected area with generous amounts of water. In case of contact with the eyes, immediately rinse with plenty of water and seek medical attention.
6. DeTrey Conditioner 36 gel should extrude easily: **DO NOT USE EXCESSIVE FORCE.** Replace original cap of DeTrey Conditioner 36 tightly after each use to avoid evaporation. Discard needle after use, as needles may clog if gel is allowed to dry inside.

## PRECAUTIONS

Avoid Prime&Bond NT saturating gingival retraction cord. If Prime&Bond NT soaks into the cord, it may set hard and bond the cord to the underlying tooth surface making removal difficult.

## INTERACTIONS WITH DENTAL MATERIALS

### DeTrey Conditioner 36:

Some liners and bases may be etched by phosphoric acid. In general, this does not impair their barrier function.

### Prime&Bond NT:

If H<sub>2</sub>O<sub>2</sub> has been used to clean the cavity, proper rinsing is essential. Higher concentration H<sub>2</sub>O<sub>2</sub> may interfere with the setting of polymerisable material and should not be used prior to the application of Prime&Bond NT.

Prolonged and intensive contact with acetone-containing products may lead to minute dissolution of the outermost surface of calcium hydroxide materials. This has no detrimental effect on the adhesion to the cavity walls.

### Prime&Bond NT/SureFil:

Eugenol containing dental materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.

## ADVERSE REACTIONS

### Prime&Bond NT:

The following adverse reaction has been associated with the use of acetone solutions and acrylate monomers:

- Reversible inflammatory changes of the oral mucosa after accidental contact.

**SureFil:**

Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.

**STEP-BY-STEP INSTRUCTIONS FOR USE****1. Cavity Preparation**

Cavity design requirements are essentially the same as for the preparation for a conventional composite restoration. Internal line and point angles should be rounded. No residual amalgam or other base material should be left in the internal forms of the preparation which would interfere with the Prime&Bond NT Adhesive bond to tooth or with light transmission and the hardening of the restorative. Rinse completed preparation thoroughly with air/water spray.

**2. Pulp Protection**

Prior to Acid-Etching in normal depth (routine) cavity preparations, where remaining dentin thickness is judged to be greater than 1 mm, a calcium hydroxide base is not required. In all deep preparations, in close proximity to the pulp, less than 1 mm remaining dentin thickness, place a calcium hydroxide liner such as Dycal® or Prisma® VLC Dycal Base/Liner Composition.

**3. Placement of Matrix**

The use of a deadsoft, thin matrix band and subsequent burnishing of the matrix band will improve final interproximal contact and contour.

PRE-WEDGING IS ADVOCATED TO ACHIEVE SLIGHT SEPARATION AND FACILITATE ACCEPTABLE PROXIMAL CONTACT.

**4. Acid Conditioning of Enamel and Dentine (Total Etch Technique)<sup>1</sup>**

When used as a bonding agent for composite materials, it is recommended to follow the Total Etch Technique described below:

**4.1 Application of DeTrey Conditioner 36**

Gently extrude DeTrey Conditioner 36 gel (36% phosphoric acid) to the cavity surfaces starting at the enamel margins. For best results, condition enamel for at least 15 seconds and dentine for 15 seconds or less.

**4.2 Rinsing and Drying**

Remove gel with aspirator tube and/or vigorous water spray and rinse conditioned areas thoroughly for at least 15 seconds.

Remove water from the rinsed cavity with a soft blow of air. Avoid desiccating the dentine, leave a moist surface.

Once the surfaces have been properly treated, they must be kept uncontaminated. If salivary contamination occurs, thoroughly clean with vigorous water-spray, dry, and repeat conditioning procedure of enamel for 5 seconds only. Rinse and dry as described above.

**5. Application of Prime&Bond NT**

One layer of Prime&Bond NT is applied:

5.1 Dispense Prime&Bond NT directly onto a fresh Applicator Tip<sup>2</sup> or onto a disposable brush. Alternatively, dispense into a fresh DENTSPLY Applicator Dish<sup>2</sup> or standard dappen dish.

<sup>1</sup> Alternatively to the Total Etch Technique, the conventional Enamel Etch Technique can be followed. In this case, the enamel margins only are treated with DeTrey Conditioner 36 for at least 15 seconds. Then rinse and dry as described below.

<sup>2</sup> DENTSPLY Applicator Dish and Applicator Tips are available from your dental dealer.

- 5.2 Immediately apply ample amounts of Prime&Bond NT to thoroughly wet all tooth surfaces. These surfaces should be saturated which may necessitate additional application of Prime&Bond NT.
- 5.3 Leave the surface undisturbed for 20 seconds.
- 5.4 Remove solvent by blowing gently with air from a dental syringe for at least 5 seconds. Surface should have a uniform, glossy appearance. If not, repeat steps 2 to 4.
- 5.5 Light-cure for a minimum of 10 seconds<sup>3</sup>. Ensure uniform exposure of all cavity surfaces.
- 5.6 Immediately place SureFil restorative over the cured Prime&Bond NT.

## **6. Placement of SureFil High Density Posterior Restorative**

- 6.1 Dispense the necessary amount of SureFil restorative material from the syringe onto a mixing pad by turning the handle slowly in a clockwise direction. To prevent oozing of the material when dispensing is completed, point the front tip of the syringe upwards and turn the handle anti-clockwise. Immediately re-close the syringe with the respective cap. Place SureFil in increments into the cavity and protect remaining material against light. SureFil may be placed up to 5 mm, light-curing for 40 seconds.
- 6.2 Carry and condense enough material to slightly (0.1 mm) overfill the cavity. Contour and shape with the operator's choice of clean amalgam or composite carving and burnishing instruments. To minimize finishing time and reduce the likelihood of ditching margins, carve margins and anatomy to final form. SureFil resists slumping, allowing carving of the majority of anatomical form prior to Visible Light Curing.

### Important Technique Note:

It is strongly recommended that the instrument used to contour cavo-surface margins (occlusal and proximal) should be lubricated from time to time with a thin coat of residual Prime&Bond NT to ensure optimal marginal adaptation. (Prime&Bond NT should be previously dispensed to allow for solvent evaporation. If using product immediately dispensed, lightly apply air to material in the plastic dappen dish for evaporation of the acetone solvent prior to use as a „modeling“ liquid.) Suggested instruments for carving may be either blunt round instruments such as PKT-3 or a round "ice cream cone type" instrument; alternative carving instruments include a cleoid discoid, Walls or Hollenback carver. Be sure the instrument is well lubricated with residual Prime&Bond NT during the carving and contour process.

- 6.3 Expose each area of the restoration surface to the Spectrum® or Prolite® Curing Unit for at least 40 seconds. The composite should be additionally exposed to the curing unit through the proximal, lingual and buccal enamel walls following metal matrix removal.

## **FINISHING AND POLISHING**

Begin finishing immediately after matrix removal and final increment curing. Gross excess may be removed with Prisma® Finishing Burs or other carbide finishing burs. Additional finishing and polishing is obtained by use of Enhance™ Discs, Cups or Points and interproximal strips. Alternatively, other standard aluminum oxide disc series and/or abrasive impregnated points may be used. A high luster can be obtained on SureFil material by applying Prisma® Gloss™ followed by Prisma Gloss Extrafine with Enhance Polishing Foam Cups.

---

<sup>3</sup> When using a high-performance unit such as the curing lights manufactured by DENTSPLY, a curing time of 10 seconds is sufficient. For curing lights with an output lower than 300 mW/cm<sup>2</sup>, a curing time of 20 seconds should be observed.

## Procedure

**STEP 1:** Complete gross reduction of excess and general outline form of the restoration using Prisma Finishing Burs, or diamond finishing instruments.

**STEP 2:** Insert an Enhance Finishing Disc, Cup or Point into a latch conventional speed contra-angled handpiece and continue finishing. The aggressiveness of the Enhance Disc is controlled by the pressure applied to the surface at the composite. The greater the pressure, the more material is removed; lighter pressure leaves a smooth surface without removing bulk.

**STEP 3:** Attach an Enhance Polishing Cup to the provided mandrel by inserting the bayonet end of the mandrel into the opening on the narrow end of the cup. Rotating the mandrel 1/4 turn makes insertion easier. Make sure the mandrel inserts fully into the cup. Insert the mandrel into a conventional speed contra-angled handpiece.

**STEP 4:** Apply a small amount of Prisma Gloss Material to the surface of the Enhance Cup. Work the surface of the restoration, dry initially, at moderate speed and pressure. Use flat end and corner edge of cup.

**STEP 5:** To increase surface luster, ADD WATER IN SMALL AMOUNTS (i.e. dropwise) to dilute the paste using a light circular buffing motion. Repeat as needed to produce a smooth surface.

**STEP 6:** Rinse Prisma Gloss from the tooth surface and Enhance Cup. Apply Prisma Gloss Extrafine Paste to the Enhance Cup as described in step 4 and polish the surfaces dry at first, then adding increasing amounts of water for 15 - 30 seconds for final luster. The Enhance Cup should be discarded after use.

### Technique Hints:

1. The aggressiveness of the Enhance Disc is controlled by the pressure applied to the surface of the composite. The greater the pressure, the more material is removed; lighter pressure leaves a smooth surface without removing bulk.
2. Polishing efficiency is greatly increased and splatter of Prisma Gloss is greatly reduced by rubbing paste into the surface of the Enhance Polishing Cup before use. This procedure leaves a thin, adherent surface film of paste in the cup surface. Use the paste impregnated cup dry at low speed; after 30 - 60 seconds, apply a small amount of water to the tooth and cup. Use again at low speed for another 30 - 60 seconds to achieve high surface luster.
3. As with any rotary instrument, heat will build up when using the Enhance Disc with prolonged contact. Use with intermittent pressure.
4. Use the flat end and corner edge of the Enhance Polishing Cup. Excessive lateral pressure may dislodge cup from mandrel.

## STORAGE

### **All products:**

Keep out of sunlight. Not to be stored at temperatures exceeding 24 °C.

### **DeTrey Conditioner 36:**

Replace cap immediately after use.

### **Prime&Bond NT:**

The Prime&Bond NT bottle should be tightly closed immediately after use. Keep in a well ventilated place.

Protect from exposure to ambient light.

### **SureFil:**

Do not freeze. Replace cap upon extrusion of composite. Keep in a well ventilated place. Protect from exposure to ambient light.

**BATCH NUMBER AND EXPIRY**

The batch number should be quoted in all correspondence which requires identification of the product.

Do not use after expiry date.

If you have any questions, please contact:

**EU-Authorized Representative:****DENTSPLY DeTrey GmbH****De-Trey-Str. 1****78467 Konstanz****GERMANY****Phone +49 (0) 75 31 5 83-0****Distributor:****DENTSPLY Limited****Hamm Moor Lane****Addlestone, Weybridge****Surrey KT15 2SE****Phone (0 19 32) 85 34 22**

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© DENTSPLY DeTrey 1998-05-14